

**Anti-TMEM87A antibody (STJA0008840)**

STJA0008840

**GENERAL INFORMATION**

|                          |                                                                                                                      |
|--------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Transmembrane protein 87A is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB                                                                                                                   |
| <b>Host/Source</b>       | Rabbit                                                                                                               |
| <b>Reactivity</b>        | Human/Mouse                                                                                                          |

**PRODUCT PROPERTIES**

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                            |
| <b>Clone ID</b>            |                                                                                                                       |
| <b>Concentration</b>       | 1 mg/mL                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:500-2000                                                                                                         |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

**TARGET INFORMATION**

|                           |                                                                    |
|---------------------------|--------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">25963</a>                                              |
| <b>Gene Symbol</b>        | <a href="#">TMEM87A</a>                                            |
| <b>Uniprot ID</b>         | <a href="#">TM87A_HUMAN</a>                                        |
| <b>Immunogen</b>          | Synthesized peptide derived from human TMEM87A                     |
| <b>Immunogen Region</b>   |                                                                    |
| <b>Specificity</b>        | This antibody detects endogenous levels of TMEM87A at Human, Mouse |
| <b>Immunogen Sequence</b> |                                                                    |