

Anti-ANKRD50 antibody (STJA0008901)

STJA0008901

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ankyrin repeat domain-containing protein 50 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57182
Gene Symbol	ANKRD50
Uniprot ID	ANR50_HUMAN
Immunogen	Synthesized peptide derived from human ANKRD50
Immunogen Region	
Specificity	This antibody detects endogenous levels of ANKRD50 at Human
Immunogen Sequence	