

Anti-PPP1R9A antibody (STJA0009076)

STJA0009076

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Neurabin-1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55607
Gene Symbol	PPP1R9A
Uniprot ID	NEB1_HUMAN
Immunogen	Synthesized peptide derived from human PPP1R9A
Immunogen Region	
Specificity	This antibody detects endogenous levels of PPP1R9A at Human, Rat
Immunogen Sequence	