

Anti-ART5 antibody (STJA0009120)

STJA0009120

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ecto-ADP-ribosyltransferase 5transferase 5+---arginine ADP-ribosyltransferase 5 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116969
Gene Symbol	ART5
Uniprot ID	NAR5_HUMAN
Immunogen	Synthesized peptide derived from human ART5
Immunogen Region	
Specificity	This antibody detects endogenous levels of ART5 at Human, Mouse, Rat
Immunogen Sequence	