

Anti-TYW1B antibody (STJA0009174)

STJA0009174

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-S-adenosyl-L-methionine-dependent tRNA 4-demethylwyosine synthase TYW1B is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	441250
Gene Symbol	TYW1B
Uniprot ID	TYW1B_HUMAN
Immunogen	Synthesized peptide derived from human TYW1B
Immunogen Region	
Specificity	This antibody detects endogenous levels of TYW1B at Human
Immunogen Sequence	