

Anti-PRRG1 antibody (STJA0009268)

STJA0009268

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Transmembrane gamma-carboxyglutamic acid protein 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5638
Gene Symbol	PRRG1
Uniprot ID	TMG1_HUMAN
Immunogen	Synthesized peptide derived from human PRRG1
Immunogen Region	
Specificity	This antibody detects endogenous levels of PRRG1 at Human
Immunogen Sequence	