

Anti-NLGN4Y antibody (STJA0009352)

STJA0009352

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Neurologin-4, Y-linked is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22829
Gene Symbol	NLGN4Y
Uniprot ID	NLGNY_HUMAN
Immunogen	Synthesized peptide derived from human NLGN4Y
Immunogen Region	
Specificity	This antibody detects endogenous levels of NLGN4Y at Human
Immunogen Sequence	