

Anti-HLCS antibody (STJA0009414)

STJA0009414

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Biotin--protein ligase ligase; Biotin--propionyl-CoA-carboxylase ligase; Biotin--methylcrotonoyl-CoA-carboxylase ligase; Biotin--acetyl-CoA-carboxylase ligase is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3141
Gene Symbol	HLCS
Uniprot ID	BPL1_HUMAN
Immunogen	Synthesized peptide derived from human HLCS
Immunogen Region	
Specificity	This antibody detects endogenous levels of HLCS at Human, Mouse
Immunogen Sequence	