

Anti-B3GNT7 antibody (STJA0009439)

STJA0009439

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UDP-GlcNAc:betaGal beta-1, 3-N-acetylglucosaminyltransferase 7 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	93010
Gene Symbol	B3GNT7
Uniprot ID	B3GN7_HUMAN
Immunogen	Synthesized peptide derived from human B3GNT7
Immunogen Region	
Specificity	This antibody detects endogenous levels of B3GNT7 at Human, Mouse, Rat
Immunogen Sequence	