

Anti-DMRT1 antibody (330-373) (STJA0009517)

STJA0009517

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Doublesex-and mab-3-related transcription factor 1 (330-373) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1761
Gene Symbol	DMRT1
Uniprot ID	DMRT1_HUMAN
Immunogen	Synthesized peptide derived from human DMRT1. AA range:330-373
Immunogen Region	330-373
Specificity	This antibody detects endogenous levels of DMRT1 at Human, Mouse
Immunogen Sequence	