

Anti-C2CD2L antibody (STJA0009539)

STJA0009539

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospholipid transfer protein C2CD2L is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9854
Gene Symbol	C2CD2L
Uniprot ID	C2C2L_HUMAN
Immunogen	Synthesized peptide derived from human TMEM24
Immunogen Region	
Specificity	This antibody detects endogenous levels of TMEM24 at Human, Mouse
Immunogen Sequence	