

Anti-ADAP1 antibody (STJA0009551)

STJA0009551

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Arf-GAP with dual PH domain-containing protein 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11033
Gene Symbol	ADAP1
Uniprot ID	ADAP1_HUMAN
Immunogen	Synthesized peptide derived from human ADAP1
Immunogen Region	
Specificity	This antibody detects endogenous levels of ADAP1 at Human
Immunogen Sequence	