

Anti-PDS5A antibody (STJA0009683)

STJA0009683

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Sister chromatid cohesion protein PDS5 homolog A is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23244
Gene Symbol	PDS5A
Uniprot ID	PDS5A_HUMAN
Immunogen	Synthesized peptide derived from human PDS5A
Immunogen Region	
Specificity	This antibody detects endogenous levels of PDS5A at Human, Rat
Immunogen Sequence	