

**Anti-Phospho-BRSK2-Thr260 antibody (STJA0009843)**

STJA0009843

**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Primary antibodies  |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Phospho-Serine/threonine-protein kinase BRSK2-Thr260 is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB/ELISA  |
| <b>Host/Source</b>       | Rabbit  |
| <b>Reactivity</b>        | Human/Mouse/Rat   |

**PRODUCT PROPERTIES**

|                       |   |
|-----------------------|---|
| <b>Clonality</b>      | Polyclonal  |
| <b>Clone ID</b>       |   |
| <b>Concentration</b>  | 1 mg/mL   |
| <b>Conjugation</b>    | Unconjugated  |
| <b>Purification</b>   | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| <b>Dilution Range</b> | WB 1:1000-2000<br>ELISA 1:5000-20000  |
| <b>Formulation</b>    | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                   |
| <b>Isotype</b>        | IgG   |
| <b>Storage</b>        | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.            |
| <b>Instruction</b>    |   |

**TARGET INFORMATION**

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | 9024  |
| <b>Gene Symbol</b>        | BRSK2   |
| <b>Uniprot ID</b>         | BRSK2_HUMAN   |
| <b>Immunogen</b>          | Synthesized peptide derived from human BRSK2 (Phospho Thr260)                       |
| <b>Immunogen Region</b>   |   |
| <b>Specificity</b>        | This antibody detects endogenous levels of Human, Mouse, Rat BRSK2 (Phospho Thr260) |
| <b>Immunogen Sequence</b> |   |