

Anti-Phospho-TYK2-Tyr292 antibody (STJA0009867)

STJA0009867

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Non-receptor tyrosine-protein kinase TYK2-Tyr292 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7297
Gene Symbol	TYK2
Uniprot ID	TYK2_HUMAN
Immunogen	Synthesized peptide derived from human Tyk 2 (Phospho Tyr292)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human Tyk 2 (Phospho Tyr292)
Immunogen Sequence	