

Anti-Phospho-TLR3-Tyr759 antibody (STJA0009953)

STJA0009953

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Toll-like receptor 3-Tyr759 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7098
Gene Symbol	TLR3
Uniprot ID	TLR3_HUMAN
Immunogen	Synthesized peptide derived from human TLR3 (Phospho-Tyr759)
Immunogen Region	
Specificity	This antibody detects endogenous levels of TLR3 (Phospho-Tyr759) at Human, Mouse, Rat
Immunogen Sequence	