

Anti-Phospho-SMAD1-Ser463+Ser465 antibody (STJA0009975)

STJA0009975

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Smad1/5/9-Ser463+Ser465 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4090 4093 4086 SMAD5 SMAD9 SMAD5_HUMAN SMAD9_HUMAN SMAD1_HUMAN
Immunogen	Synthesized peptide derived from human Smad1/5/9 (Phospho-Ser463+Ser465)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Smad1/5/9 (Phospho-Ser463+Ser465) at Human, Mouse, Rat
Immunogen Sequence	



Western Blot analysis of various, using primary antibody at 1:1000 dilution. Secondary antibody (STJS000791) was diluted at 1:10000