

Anti-Phospho-NACA-Thr159 antibody (STJA0010000)

STJA0010000

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Nascent polypeptide-associated complex subunit alpha-Thr159 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4666
Gene Symbol	NACA
Uniprot ID	NACA_HUMAN
Immunogen	Synthesized peptide derived from human NACA (Phospho-Thr159)
Immunogen Region	
Specificity	This antibody detects endogenous levels of NACA (Phospho-Thr159) at Human, Mouse, Rat
Immunogen Sequence	