

Anti-ACAP1 antibody (490-570) (STJA0010088)

STJA0010088

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Arf-GAP with coiled-coil, ANK repeat and PH domain-containing protein 1 (490-570) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9744
Gene Symbol	ACAP1
Uniprot ID	ACAP1_HUMAN
Immunogen	Synthesized peptide derived from human Centaurin-Beta 1 AA range: 490-570
Immunogen Region	490-570
Specificity	This antibody detects endogenous levels of Human Centaurin-Beta 1
Immunogen Sequence	