

Anti-HMGN1 antibody (1-50) (STJA0010135)

STJA0010135

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-HMGN1/2/3/4 (1-50) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9324 3151 10473 HMGN3 HMGN2 HMGN3_HUMAN HMGN2_HUMAN HMGN4_HUMAN
Immunogen	Synthesized peptide derived from human HMGN1/2/3/4 AA range: 1-50
Immunogen Region	1-50
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat HMGN1/2/3/4
Immunogen Sequence	