

Anti-IL17A antibody (50-100) (STJA0010140)

STJA0010140

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Interleukin-17A (50-100) is suitable for use in Western Blot, Flow Cytometry and ELISA research applications.
Applications	WB/FC/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 Flow Cyt 1:50-200 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3605
Gene Symbol	IL17A
Uniprot ID	IL17_HUMAN
Immunogen	Synthesized peptide derived from human IL-17 AA range: 50-100
Immunogen Region	50-100
Specificity	This antibody detects endogenous levels of Human IL-17
Immunogen Sequence	