

Anti-CNN1 antibody (111-160) (STJA0010156)

STJA0010156

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Calponin-1 (111-160) is suitable for use in Western Blot and ELISA research applications. |
| Applications | WB/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Rat/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:1000-2000 ELISA 1:5000-20000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 1264 |
| Gene Symbol | CNN1 |
| Uniprot ID | CNN1_HUMAN |
| Immunogen | Synthesized peptide derived from human Calponin 1 AA range: 111-160 |
| Immunogen Region | 111-160 |
| Specificity | This antibody detects endogenous levels of Human Calponin 1 |
| Immunogen Sequence | |