

Anti-SDCBP antibody (71-120) (STJA0010163)

STJA0010163

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Syntenin-1 (71-120) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Dog/Pig

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6386
Gene Symbol	SDCBP
Uniprot ID	SDCB1_HUMAN
Immunogen	Synthesized peptide derived from human Syntenin-2 AA range: 71-120
Immunogen Region	71-120
Specificity	This antibody detects endogenous levels of Human, Rat, Dog, Pig Syntenin-2
Immunogen Sequence	