

Anti-ATG7 antibody (521-570) (STJA0010166)

STJA0010166

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ubiquitin-like modifier-activating enzyme ATG7 (521-570) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10533
Gene Symbol	ATG7
Uniprot ID	ATG7_HUMAN
Immunogen	Synthesized peptide derived from human APG7 AA range: 521-570
Immunogen Region	521-570
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat APG7
Immunogen Sequence	