

Anti-DLG4 antibody (530-610) (STJA0010242)

STJA0010242

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DLG4 (530-610) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1742
Gene Symbol	DLG4
Uniprot ID	DLG4_HUMAN
Immunogen	Synthesized peptide derived from human DLG4 AA range: 530-610
Immunogen Region	530-610
Specificity	This antibody detects endogenous levels of Human, Mouse, Rat DLG4
Immunogen Sequence	