

Anti-PPM1G antibody (C-Term) (STJA0010311)

STJA0010311

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein phosphatase 1G (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5496
Gene Symbol	PPM1G
Uniprot ID	PPM1G_HUMAN
Immunogen	Synthesized peptide derived from human C-terminal PPM1G
Immunogen Region	C-Term
Specificity	This antibody detects endogenous levels of PPM1G at Human, Mouse
Immunogen Sequence	