

Anti-PGK1 antibody (N-Term) (STJA0010315)

STJA0010315

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phosphoglycerate kinase 1 (N-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5230
Gene Symbol	PGK1
Uniprot ID	PGK1_HUMAN
Immunogen	Synthesized peptide derived from human N-terminal PGK1
Immunogen Region	N-Term
Specificity	This antibody detects endogenous levels of PGK1 at Human, Mouse, Rat
Immunogen Sequence	