

Anti-HMGCL antibody (C-Term) (STJA0010324)

STJA0010324

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Hydroxymethylglutaryl-CoA lyase, mitochondrial (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3155
Gene Symbol	HMGCL
Uniprot ID	HMGCL_HUMAN
Immunogen	Synthesized peptide derived from human C-terminal HMGCL
Immunogen Region	C-Term
Specificity	This antibody detects endogenous levels of HMGCL at Human, Mouse, Rat
Immunogen Sequence	