

Anti-CUL9 antibody (C-Term) (STJA0010340)

STJA0010340

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cullin-9 (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23113
Gene Symbol	CUL9
Uniprot ID	CUL9_HUMAN
Immunogen	Synthesized peptide derived from human C-terminal CUL9
Immunogen Region	C-Term
Specificity	This antibody detects endogenous levels of CUL9 at Human, Mouse
Immunogen Sequence	