

Anti-WWP1 antibody (N-Term) (STJA0010355)

STJA0010355

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-NEDD4-like E3 ubiquitin-protein ligase WWP1 (N-Term) is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC/WB |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:500-2000 IHC 1:50-200 |
| Formulation | Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 11059 |
| Gene Symbol | WWP1 |
| Uniprot ID | WWP1_HUMAN |
| Immunogen | Synthesized peptide derived from human N-terminal WWP1 |
| Immunogen Region | N-Term |
| Specificity | This antibody detects endogenous levels of WWP1 at Human, Mouse |
| Immunogen Sequence | |