

Anti-MMP9 antibody [SM0375] (STJA0010375)

STJA0010375

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-MMP9 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0375
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H) 1:300-1:600
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG2ak
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4318
Gene Symbol	MMP9
Uniprot ID	MMP9_HUMAN
Immunogen	Recombinant protein corresponding to Human MMP9
Immunogen Region	
Specificity	
Immunogen Sequence	