

Anti-Phospho-JNK-T183/Y185 antibody [SM0384] (STJA0010384)

STJA0010384

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Phospho-JNK-T183/Y185 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC/IF
Host/Source	Mouse
Reactivity	Human





PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0384
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H) 1:500-1:1000 IHC/IF (H) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG2bk
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID [5599](#)
[5602](#)
[5601](#)
[MAPK8](#)
[MAPK10](#)
[MK08_HUMAN](#)
[MK10_HUMAN](#)
[MK09_HUMAN](#)

Immunogen KLH conjugated Synthetic peptide corresponding to Human JNK (T183/Y185)
Immunogen Region
Specificity
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081