

Anti-V5-Tag antibody [SM0385] (STJA0010385)

STJA0010385

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-V5-Tag is suitable for use in Western Blot, Immunocytochemistry and Immunofluorescence research applications.
Applications	WB/ICC/IF
Host/Source	Mouse
Reactivity	V5-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0385
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	KLH conjugated Synthetic peptide corresponding to V5-tag
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081