

Anti-NAPSIN A antibody [SM0409] (STJA0010409)

STJA0010409

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-NAPSIN A is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0409
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M, R) 1:500-1:1500
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG2b
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9476
Gene Symbol	NAPSA
Uniprot ID	NAPSA_HUMAN
Immunogen	Recombinant protein corresponding to Human NAPSIN-A
Immunogen Region	
Specificity	
Immunogen Sequence	