

Anti-IL-6 antibody [SM0423] (STJA0010423)

STJA0010423

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-IL-6 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0423
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M) 1:300-1:600
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16193
Gene Symbol	IL6
Uniprot ID	IL6_MOUSE
Immunogen	Recombinant protein corresponding to Mouse IL-6
Immunogen Region	
Specificity	
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081