

Anti-Peroxiredoxin 2/PRP antibody [SM0428] (STJA0010428)

STJA0010428

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Peroxiredoxin 2/PRP is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Mouse |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | SM0428 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB (H, M) 1:500-1:1000 |
| Formulation | PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol. |
| Isotype | IgG |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 21672 |
| Gene Symbol | Prdx2 |
| Uniprot ID | PRDX2_MOUSE |
| Immunogen | KLH conjugated Synthetic peptide corresponding to Mouse PRP |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081