

## Anti-Cyclophilin A antibody [SM0429] (STJA0010429)

STJA0010429

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Cyclophilin A is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SM0429
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB (H, M, R) 1:500-1:1000
<b>Formulation</b>	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
<b>Isotype</b>	IgG1k
<b>Storage Instruction</b>	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">268373</a>
<b>Gene Symbol</b>	<a href="#">Ppia</a>
<b>Uniprot ID</b>	<a href="#">PPIA_MOUSE</a>
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Mouse Cyclophilin A
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	