

Anti-CDKN2A/p16INK4a antibody [SM0438] (STJA0010438)

STJA0010438

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CDKN2A/p16INK4a is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0438
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H) 1:500-1:2000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1029
Gene Symbol	CDKN2A
Uniprot ID	CDN2A_HUMAN
Immunogen	KLH conjugated Synthetic peptide corresponding to Human P16
Immunogen Region	
Specificity	
Immunogen Sequence	