

**Anti-iNOS antibody [SM0487] (STJA0010487)**

STJA0010487

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-iNOS is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SM0487
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB (M) 1:500-1:1000
<b>Formulation</b>	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
<b>Isotype</b>	IgG2bk
<b>Storage Instruction</b>	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">4843</a>
<b>Gene Symbol</b>	<a href="#">NOS2</a>
<b>Uniprot ID</b>	<a href="#">NOS2_HUMAN</a>
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Mouse iNOS
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	

---

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081