

**Anti-Recombinant-Androgen receptor antibody [SR0520] (STJA0010520)**

STJA0010520

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit monoclonal antibody anti-Recombinant-Androgen receptor is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	IHC/IF
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SR0520
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Recombinant Expression and Affinity purified
<b>Dilution Range</b>	IHC/IF (H, M, R) 1:200-1:400
<b>Formulation</b>	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
<b>Isotype</b>	
<b>Storage</b>	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">367</a>
<b>Gene Symbol</b>	<a href="#">AR</a>
<b>Uniprot ID</b>	<a href="#">ANDR_HUMAN</a>
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Human Androgen Receptor
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	