

Anti-Recombinant-SOX9 antibody [SR0521] (STJA0010521)

STJA0010521

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-SOX9 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0521
Concentration	
Conjugation	Unconjugated
Purification	Recombinant Expression and Affinity purified
Dilution Range	IHC/IF (H, M, R) 1:200-1:400
Formulation	PBS, 50% glycerol, 0.05% Proclin 300, 0.05%BSA
Isotype	
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6662
Gene Symbol	SOX9
Uniprot ID	SOX9_HUMAN
Immunogen	KLH conjugated Synthetic peptide corresponding to Human SOX9
Immunogen Region	
Specificity	
Immunogen Sequence	