

Anti-Recombinant-Phospho-STAT3-Y705 antibody [SR0523] (STJA0010523)

STJA0010523

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit monoclonal antibody anti-Recombinant-Phospho-STAT3-Y705 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications. |
| Applications | WB/IHC/IF |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse/Rat |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | SR0523 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | WB (H, M, R) 1:1000-1:2000 IHC/IF (H, M, R) 1:500-1:2000 |
| Formulation | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol. |
| Isotype | IgG |
| Storage | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 6774 |
| Gene Symbol | STAT3 |
| Uniprot ID | STAT3_HUMAN |
| Immunogen | KLH conjugated Synthetic phosphopeptide corresponding to Human STAT3 (Tyr705) |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |