

Anti-Recombinant-CD34 antibody [SR0541] (STJA0010541)

STJA0010541

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-CD34 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0541
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M, R, Mk) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	947
Gene Symbol	CD34
Uniprot ID	CD34_HUMAN
Immunogen	KLH conjugated Synthetic peptide corresponding to Human CD34
Immunogen Region	
Specificity	
Immunogen Sequence	