

Anti-Recombinant-Phospho-CA125/MUC16 antibody [SM0545] (STJA0010545)

STJA0010545

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Phospho-CA125/MUC16 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0545
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H) 1:250-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	94025
Gene Symbol	MUC16
Uniprot ID	MUC16_HUMAN
Immunogen	Human recombinant CA125 protein expressed in HEK293 cells
Immunogen Region	
Specificity	
Immunogen Sequence	