

Anti-Recombinant-Phospho-V5-Tag antibody [SM0547] (STJA0010547)

STJA0010547

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Phospho-V5-Tag is suitable for use in Western Blot, Immunocytochemistry and Immunofluorescence research applications.
Applications	WB/ICC/IF
Host/Source	Mouse
Reactivity	V5-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0547
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1000 ICC/IF (Species independent) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	KLH conjugated Synthetic peptide corresponding to V5 tag
Immunogen Region	
Specificity	
Immunogen Sequence	