

Anti-Recombinant-Desmin antibody [SM0557] (STJA0010557)

STJA0010557

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-Recombinant-Desmin is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications. |
| Applications | WB/IHC/IF |
| Host/Source | Mouse |
| Reactivity | Mouse/Rat |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Monoclonal |
| Clone ID | SM0557 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC/IF (M, R) 1:500-1:1000 |
| Formulation | PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol. |
| Isotype | IgG1 |
| Storage | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |
| Instruction | |

TARGET INFORMATION

| | |
|--------------------|---|
| Gene ID | 13346 |
| Gene Symbol | Des |
| Uniprot ID | DESM_MOUSE |
| Immunogen | Recombinant protein corresponding to Mouse desmin |
| Immunogen | |
| Region | |
| Specificity | |
| Immunogen | |
| Sequence | |