

Anti-Recombinant-Granzyme B antibody [SM0563] (STJA0010563)

STJA0010563

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Granzyme B is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0563
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3002
Gene Symbol	GZMB
Uniprot ID	GRAB_HUMAN
Immunogen	Human recombinant Granzyme B protein expressed in HEK293 cells
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	