

## Anti-Recombinant-GAP43 antibody [SR0564] (STJA0010564)

STJA0010564

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit monoclonal antibody anti-Recombinant-GAP43 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	WB/IHC/IF
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SR0564
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB (M, R) 1:500-1:1000 IHC/IF (M, R) 1:500-1:1000
<b>Formulation</b>	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">14432</a>
<b>Gene Symbol</b>	<a href="#">Gap43</a>
<b>Uniprot ID</b>	<a href="#">NEUM_MOUSE</a>
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Mouse GAP43
<b>Immunogen</b>	
<b>Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081