

Anti-Recombinant-Mast cell tryptase antibody [SM0576] (STJA0010576)

STJA0010576

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Mast cell tryptase is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0576
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7177
Gene Symbol	TPSAB1
Uniprot ID	TRYB1_HUMAN
Immunogen	Recombinant protein corresponding to Human Mast cell tryptase
Immunogen Region	
Specificity	
Immunogen Sequence	