

Anti-Recombinant-non-muscle Myosin IIB/MYH10 antibody [SM0587] (STJA0010587)

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-non-muscle Myosin IIB/MYH10 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC/IF
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0587
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000 IHC/IF (H, M, R) 1:250-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1k
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	77579
Gene Symbol	Myh10
Uniprot ID	MYH10_MOUSE
Immunogen	KLHconjugatedSyntheticpeptidecorrespondingtoMouseMYH10
Immunogen Region	
Specificity	
Immunogen Sequence	